

Product Name :
Vapor Liquid Equilibrium Apparatus

Product Code :
NLAB-EDUCATIONTVET16030

Description :

Vapor Liquid Equilibrium Apparatus

Technical Specification :

Vapor Liquid Equilibrium Apparatus -

We have a range of such equipment that enhances the client's trust for the requirement of high-quality product. Vapor-Liquid Equilibrium ensures the equilibrium between the liquid and the vapor.

Any substance which is at this equilibrium known as Saturated Fluid.

It is a condition where the rate of evaporation (liquid changing to vapor) equals the rate of condensation.

School Lab Instrument is Vapor-Liquid Equilibrium Apparatus manufacturer and Supplier in India as well as worldwide.

This condition is applicable for a molecular level.

Specification:-

Frequency: 50/60Hz

Patient airway pressure range: -2kPa-10kPa

Input voltage: 100-240V AC

Tidal volume range (V t): 0 mL to 1,500 mL

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003