

**Product Name :**  
Verification of Kirchoffs Law of Voltage and Current Theorem Digital

**Product Code :**  
NLAB-ELECTRONICSAB490001



**Description :**

Verification of Kirchoffs Law of Voltage and Current Theorem Digital

**Technical Specification :**

SCOPE OF LEARNING:

VERIFICATION OF KIRCHOFF'S LAW OF VOLTAGE AND CURRENT THEOREM

TECHNICAL SPECIFICATIONS:

Digital Meters:

Ammeter 50mA DC.

Voltmeter 12V DC.

Power Supplies:

DC Supply IC Regulated 0-12V DC, 150mA.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Resistances 2Watts

Voltage Control through Potentiometer.

**SALIENT FEATURES:**

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 4mm Colored Sockets.

Patch Cords 4mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003