

Product Name :
0°C Walkin Cold Room

Product Code :
NLAB-SCIENTIFICEQ35002



Description :

0°C Walkin Cold Room

Technical Specification :

The 0°C walk-in cold rooms from Weiber are designed to provide and sustain working temperatures as low as 0°C for the preservation and storage of a wide variety of storage items including poultry, meat, fruits and vegetables, flowers, bakery items, medicines, chemicals etc. These cooling equipments provide a positive working temperature range for the bulk storage of fruits and vegetables, vaccines etc. and a negative working temperature range for the storage of meat products and other temperature sensitive materials. Weiber 0°C walk-in cold room Manufacturers provide these chambers in various sizes and configurations and can be custom-designed to suit specific requirements. Owing to their large storage capacities and flexible design, these equipments find widespread usage in butcheries, bakeries, restaurants, research laboratories, pharmacies, hospitals, mortuaries, flower shops etc.

Weiber 0 Degree Cold Room are designed for excellent insulation and refrigeration capabilities. They are made up of prefabricated PUF insulated panels of varying thickness and use pre-coated GI sheets or stainless steel sheets for external and internal surfaces. The stainless steel construction adds rigidity and structure to the equipment and ensures non-corrosive, smooth and uniform external and internal surfaces. The excellent insulation panel provides maximum thermal efficiency and reduces the power consumption. All the walls, ceiling and door are well insulated for minimum energy consumption. The Weiber Cold Room 0 Degree C are available as easy to install, light weight equipments that can be readily assembled using cam-locking arrangement.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003