

**Product Name :**  
1 Channel ECG Machine

**Product Code :**  
NLABMEDICAL-M3380005

**Description :**

1 Channel ECG Machine

**Technical Specification :**

High-resolution thermal printing array system.

Frequency response is as high as 150Hz.

Capable of printing continuously on channel trace and annotations including lead mark and parameter such as sensitivity, paper speed and filter operation status.

High resolution digital filter inhibits baseline drift without affecting ECG waveform.

One-touch operation.

Concept of floating input circuit meets IEC safety standard.

Rechargeable battery supporting more than 50 patient examinations lead Standard: 12 leads.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003