

| | |
|--|---|
| Product Name : 100MA Mobile X Ray Machines | Product Code : NLABMEDICAL-M3350003 |
|--|---|

| |
|--|
| |
|--|

Description :

100MA Mobile X Ray Machines

Technical Specification :

100MA Mobile X Ray Machines

Counter-balanced. tube arm for increased reach better positioning.

Integrated cassette box for easier operation and better cassette handling.

Elegant digital display and extremely easy to use controls.

MA meter & voltmeter for line and exposure monitoring.

Automatic calculation of mass with digital display.

Digital display of KVP.

6 Led indicator for line, 0/L, X Ray, Stand by, Bucky & fuse.

Precision solid state electronic time.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003