

Product Name :
1080P Full HD Wi-Fi Digital Microscope

Product Code :
NLAB-EDUSCHOOLEQ1730103

Description :

1080P Full HD Wi-Fi Digital Microscope

Technical Specification :

Light Source: 6 LED.

Image Sensor: 3 Mega Pixels.

Focus Range: Manual focus from 0mm to 50mm.

Video Format: H.264.

Photo Format: JPEG.

Charging power: 5V DC.

Magnification Ratio: 10x ~ 220x.

Video Capture Resolution: 1920 x 1080 (1080P), 1280 x 720 (720P), 640 x 480 (VGA).

Weight: 2lbs.

Product Size: (L)6.06 in x (R)1.73 in.

iOS/Android Interface: USB 2.0.

Still Capture Resolution: 3M, 5M, 8M, 12M.

Power Source: Rechargeable Li-ion battery (Working Time: 3hrs Charging Time: 2hrs)

Frame Rate: Max 30f/s under 600 Lus Brightness.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003