

Product Name :
12 Centrifuge

Product Code :
NLAB-EDUSCHOOLEQ45027



Description :

12 Centrifuge

Technical Specification :

12 eliminates the need to change rotors when switching between microtubes and PCR strips.

With a fixed speed that produces 2,000 xg, this centrifuge is perfect for quick spin downs.

Open the lid, and the rotor quickly decelerates for removal of samples.

The included, unique Combi-Rotor is all that is required for running 12 microtubes and/or 4 PCR strips.

Its truly is a personal centrifuge with unmatched capacity and flexibility.

Simply close the lid and the unit quickly ramps up to 5500 rpm.

At just 5.5 inches wide and less than 4.5 inches tall.

5,500 rpm / 2,000 xg	
----------------------	--

Capacity 32 x 0.2ml tubes 4 x PCR Strips	
Dimensions (W x D x H) 14.6 x 20.3 x 10.8 cm	
Weight: 1 kg	
Electrical Voltage 100-240V, 50-60 Hz, 250W	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003