

**Product Name :**  
12 Centrifuge

**Product Code :**  
NLAB-EDUSCHOOLEQ45027



**Description :**

12 Centrifuge

**Technical Specification :**

12 eliminates the need to change rotors when switching between microtubes and PCR strips.

With a fixed speed that produces 2,000 xg, this centrifuge is perfect for quick spin downs.

Open the lid, and the rotor quickly decelerates for removal of samples.

The included, unique Combi-Rotor is all that is required for running 12 microtubes and/or 4 PCR strips.

Its truly is a personal centrifuge with unmatched capacity and flexibility.

Simply close the lid and the unit quickly ramps up to 5500 rpm.

At just 5.5 inches wide and less than 4.5 inches tall.

5500 rpm / 2,000 xg

Capacity	
32 x 0.2ml tubes	
4 x PCR Strips	
Dimensions (W x D x H)	
14.6 x 20.3 x 10.8 cm	
Weight	
2.1kg	
Electrical Voltage	
100-240V, 50-60 Hz, 250W	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003