

**Product Name :**  
12 Channel ECG Machine

**Product Code :**  
NLABMEDICAL-M3380004

**Description :**

12 Channel ECG Machine

**Technical Specification :**

Simultaneous acquisition of 12-lead ECG data, real-time recording and printing of 12-channel ECG waveform.

Graphic display of 12-lead ECG waveform.

Built-in analysis software of age which assures accurate analysis result.

Auto-measurement, auto-interpretation, waveform playback and storage of ECG data.

R5232 and USB interface cope with data-share or remote data management requirement.

Option wireless function makes more convenient communication with PC.

User friendly operation system contributed by built-in HELP function.

Literal and graphic operation interface.

Powerful filters to minimize interference.

Heart rate measurement and pace-maker protection circuit.

Multi-printing formats: manual/automatic, standard 12 channel, 3 channel plus 3 rhythm lead, 6 channel, 6 channel plus rhythm lead, 60s analysis of arrhythmia, R-R histogram, trend graph.

AC, DC or built-in lithium battery power supply, alarm of battery weak and lead-off.

Comfort and durability contributed by silicon gel keypad.

Tremendous ECG data can be saved in built-in SD card ECG 330G.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003