

Product Name :
12 Value Capacitor Box

Product Code :
NLAB-SCHOOLABS13002



Description :

12 Value Capacitor Box

Technical Specification :

Contains values $\times .022 \mu\text{F}$, $.01 \mu\text{F}$, $.0047 \mu\text{F}$, $.0022 \mu\text{F}$, $.001 \mu\text{F}$, 100 pF Values $\times .033 \mu\text{F}$, $.47 \mu\text{F}$, $0.1 \mu\text{F}$, $.22 \mu\text{F}$, $.33 \mu\text{F}$, $.47 \mu\text{F}$

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003