

Product Name :
12 Value Resistance Box

Product Code :
NLAB-SCHOOLABS13001



Description :

12 Value Resistance Box

Technical Specification :

Contains values ??4.7 K, 2.2 K, 1.0 K, 470 ? ,220 ? , 100 ? values ??10 K, 22 K , 47K, 100K, 220K, 470K

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003