

Product Name :
1600 MI Membrane Densitometer

Product Code :
NLAB-EDUCATIONTVET15060



Description :

1600 MI Membrane Densitometer

Technical Specification :

According to ASTM D2167

To determine the density in place of compacted or firm materials. The equipment includes a 1600 cm³ graduated column housed in an aluminum frame, a pear-type pump with shut-off valve, a base plate and a pack of 12 membranes

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003