

**Product Name :**  
3-4 Thread Overlock Machine

**Product Code :**  
NLAB-VOCATIONALTVET22001



**Description :**

3-4 Thread Overlock Machine for Textile Vocational Technology Lab

**Technical Specification :**

**Multi-Function 3-4 Thread Overlock Machine for Textile Training Lab and Textile Technology Workshop**

3 or 4 Thread, Cut and Sew Overlock Stitches

Trim seams and overcast edges

4 color threading guide

Perfect overlock stitches on all types and weights of fabric

Stitch width 5.0 - 7.0 mm

Presser foot pressure adjustment

Blind Stitch Foot for finishing and hemming in one operation

Common snap on presser foot

Toolkit

Differential feed for perfect seams on virtually any fabric type eliminates the wavy or stretched appearance of stitches on fabric

Retractable knife

Free arm / Flatbed

Trim trap for catching scraps of fabric and storage of the foot controller when the machine is not being used

Uses standard sewing machine needles

Fast and Simple Lower Looper threading system

Bright LED light to illuminate the sewing area

Thread twist protection and needle plate separation for improved sewing quality

Lightweight, compact machine with carrying handle

Built-in accessory storage

Dust cover

Instructional DVD

**Sewing:-**

Standard Presser Feet: Yes

MAX Sewing Speed: 1300spm

Stitch Length: 0.8-4.0mm

Stitch Width: 3.0-7.0mm

Foot Controller: Yes

**Display:-**

Lighting: LED

Stitch Select Display: Dial selection

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003