

**Product Name :**  
300 MHz Digital Storage Oscilloscope 2GS/s, 2 Ch, 7 Inch

**Product Code :**  
NLAB-TESTINGLABM20002



**Description :**

300 MHz Digital Storage Oscilloscope 2GS/s, 2 Ch, 7 Inch

**Technical Specification :**

Digital Storage Oscilloscope Dual Channel

Dual channel – channel bandwidth: 70MHz ~ 300MHz

The volume exquisite and it is convenient for carrying

Use Colored TFT LCD, the waveform display is clearer and stable

The real time sampling rate is 2GS/s

Memory depth is 24k Trigger Modes: Edge, Pulse, Video, Slope, Alternative and Delay trigger modes

Front and Back USB Host/Device, GPIB

Unique Digital Filter function and Waveform recorder function

Waveform recorder function

Pass/Fail function

Thirty two parameters Auto measure function

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003