

Product Name :
3000 MI Membrane Densitometer

Product Code :
NLAB-EDUCATIONTVET13003



Description :

3000 MI Membrane Densitometer

Technical Specification :

Complies with NF P94?061?2

Supplied with 6 rubber membranes

Installation accessories:

Scraper

Cutting tool

Pointed steel rod

Density spoon

Rubber mallet

Hammer 300g

Club hammer 2 kg

Peak density

Long chisel. 300 mm x width. 25 mm

Tin can with lid

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003