

Product Name :
5-Gas Analyzer

Product Code :
TK766



Description :

5-Gas Analyzer

Technical Specification :

5-Gas Analyzer

A truly portable exhaust gas analyzer that sets new standards of convenience and ease of use. It fits into one hand and can run for up to four hours on its internal re-chargeable battery.

Back lit display is four or eight lines programmable by the user.

Measurements are shown in real time.

An energy saving separate pump switch calibrates the analyzer on daily startup (run pump only during testing)

Readings are actual, not calculated.

No Special software is necessary.

This self-powered analyzer is easy to use, the only connection to the vehicle is the tailpipe.

Integral flow checking provides warning of probe blockages.

The calibration procedure uses standard automotive test gases.

Measures five gases individually:

Carbon Monoxide (CO)

Hydrocarbons(HC)

Carbon Dioxide (CO₂)

Oxygen (O₂)

Nitric Oxide (NOX)

Features:

Printer port.

Ruggedized design sets unsurpassed accuracy and durability standards for a hand held instrument.
Programmed to evaluate gasoline, LPG and CNG fuels.
Includes protective rubber boot, water trap and filter, carry case, battery charger, exhaust probe and five extra filters.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003