

**Product Name :**  
5 Liter Wet Grinder

**Product Code :**  
NLAB-VOCATIONALTVET270113

**Description :**

5 Liter Wet Grinder

**Technical Specification :**

Capacity (Litres) : 5.0 Litres, 7.0 Litres, 10.0 Litres, 15.0 Litres

Type Table Top

Material : Stainless Steel

Automation Grade : Automatic

Voltage (V) : 210 V

Frequency (Hz) : 60 Hz

Phase : 3 Phase

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003