

Product Name :
50 Kn Press for Marshall Test

Product Code :
NLAB-EDUCATIONTVET11010



Description :

50 Kn Press for Marshall Test

Technical Specification :

Compression capacity 50 kN

Speed of 50.8 mm / min

Electronic displacement sensor 50 mm stroke

Load cell 50 kN

Microprocessor-based acquisition unit with the possibility of exporting data to a computer

Supplied complete with aluminum stability mold

Marshall test software

Accessories for indirect tension / splitting test:

Splitting device for 4 "and 6" samples

Set of 2 displacement sensors

Additional displacement sensor

Software for indirect tension test

Recent configuration computer, supplied with 500VA inverter

Accessories for water sensitivity test:

Thermostatic bath with cooling group

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003