

Product Name :
6 Slice Toaster

Product Code :
NLAB-VOCATIONALTVET27017

Description :

6 Slice Toaster

Technical Specification :

6 Slice Toaster Used for making delicious & crispy toasts, this toaster is extensively used in the kitchens of homes, hotels, and restaurants.

Specifications:-

Dimension: 460*210*225 mm

Weight: 5kg

Capacity: 6 Slice

Body Material: Stainless Steel

Frequency: 60 Hz

Voltage: 210 V

Features :

Compact design

Low power consumption

Long working life

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

