

**Product Name :**  
60 MHz Digital Storage Oscilloscope 500MS s, 2 Ch, 5.7 Inch

**Product Code :**  
NLAB-TESTINGLABM20003

**Description :**

60 MHz Digital Storage Oscilloscope 500MS s, 2 Ch, 5.7 Inch

**Technical Specification :**

Dual channel – channel bandwidth: 60MHz ~ 200MHz

The volume exquisite and it is convenient for carrying

Use Colored TFT / LCD, the waveform display is clearer and stable

The real time sampling rate is 500MSa/s

Memory depth is 4K

Trigger Modes: Edge, Pulse, Video, Slope, Alternative and Delay trigger modes

Front and Back USB Host/Device, GPIB

Unique Digital Filter function and Waveform recorder function

Waveform recorder function

Pass/Fail function

Thirty two parameters Auto measure function

Cursor measure covers Manual mode, Track mode and Auto mode

Channel waveform and it's FFT waveform display on split screen

Pop-up Menu display pattern is more convenient for users using it

Multiple Language User Interface 16. Support multilingual online help

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003