

<b>Product Name :</b> 60mA 100mA X-Ray Equipments	<b>Product Code :</b> NLABMEDICAL-M3350007
<b>Description :</b>	
60mA 100mA X-Ray Equipments	
<b>Technical Specification :</b>	
Full wave rectification which delivers almost double the output than the self rectified X - Ray Systems.	
Full wave rectification with high voltage silicon rectifier for low losses and consistent radio graphic results.	
Smoothing circuit provides maximum mA output accuracy at all time settings.	
Full field light beam Collimator for proper positioning and preventing accidental cut off the field	
Continuous fluroscopic timer by pressing F_RST push button. (100mA).	
Low output power station 4Kw, for versatile, in absence of power-Line. (100mA).	
Double mA range meter with lamp. (100mA).	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003