

Product Name :
68000 Motorola Microprocessor Trainer India

Product Code :
NLAB-ENGINEERINGLB54040



Description :

68000 Motorola Microprocessor Trainer

Technical Specification :

The trainer is a valuable aid for teaching in academic institutions, software and hardware development in research institutions and R & D labs.

68000 Motorola Microprocessor Trainer is a powerful and cost-effective single-board microcomputer based on Motorola 68000 CPU supported by Motorola peripheral chips.

We are well-known manufacturers, OEM suppliers of 68000 Motorola Microprocessor Trainer for Applied Machines. Contact us for high quality 68000 Motorola Microprocessor Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003