

Product Name :
6811 Motorola Microprocessor Trainer

Product Code :
NLAB-ENGINEERINGLB54039



Description :

6811 Motorola Microprocessor Trainer

Technical Specification :

This design is ideal for experimental general experiments and also quite suited to robotics projects.

6811 is a compact trainer that is based on the 52-pin Motorola processor, and includes 32K of battery-backed static RAM, four outputs for DC motors, a connector system that allows active sensors to be individually plugged into the board, an LCD screen, and an integrated, rechargeable battery pack.

We are well-known manufacturers, OEM suppliers of 6811 Motorola Microprocessor Trainer for Applied Machines. Contact us for high quality 6811 Motorola Microprocessor Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003