

Product Name :
8 Key Differential Counter

Product Code :
NLABMEDICAL-M3100001



Description :

8 Key Differential Counter

Technical Specification :

Differential Cell Counter is commonly known as Blood Cell Counter

Blood Cell Counter is a convenient table model that can be operated with hand manually.

It is having 8 keys and 1 totalizer unit for use with a microscope to count the Basophil, Monocytes, Eosinophil, Myelofibrosis, Juven, Stab Cell, Segmental Neutrophils, Lymphocytes for other miscellaneous cells.

Each key can record up to 999 with the exception of the totalizer unit which records the total number of strokes made on all the other keys.

At totals of 100, 200, 300 and so on the totalizer unit rings a bell.

The keys are positioned low which allows the wrist to rest on the table and registers are located above the corresponding key for continuous visibility.

A reset knob at either end of the counter can turn all the keys to zero.

Blood Cell Counter has with push button

Differential Cell Counter is rust proof, water proof.

We are Manufacturer, Exporters and Suppliers of 8 Key Differential Counter or Blood Cell Counter at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003