

Product Name :
8051 Microcontroller Trainer India

Product Code :
NLAB-ENGINEERINGLB54041



Description :

8051 Microcontroller Trainer

Technical Specification :

This has created the need for a trainer and development tools.

The 8051 is a complete solution for this requirement.

8051 Microcontroller Trainer is increasing in popularity for instrumentation and control applications due to its speed and powerful instruction set – essentials for real time applications.

It is as a flexible instructional aid for academic institutions and a powerful development kit in R&D Labs.

We are well-known manufacturers, OEM suppliers of 8051 Microcontroller Trainer for Applied Machines. Contact us for high quality 8051 Microcontroller Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003