


Product Name :
8088 Microprocessor Trainer

Product Code :
NLAB-ENGINEERINGLB54038



Description :

8088 Microprocessor Trainer

Technical Specification :

The basic system can be easily expanded through the System Bus connectors.

The system is versatile and can be configured in a variety of ways via jumper options, to suit specific user requirements.

8088 Microprocessor Trainer is a powerful, general purpose micro computer trainer which can be operated either with 8086 CPU or 8088 CPU.

The CPU is configured for operating in the Max mode.

On-board provision exists for installing 8087 Numeric Data Processor.

We are well-known manufacturers, OEM suppliers of 8088 Microprocessor Trainer for Applied Machines. Contact us for high quality 8088 Microprocessor Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003