

Product Name :
A.C Ammeter 5A

Product Code :
ES251



Description :

A.C Ammeter 5A

Technical Specification :

Feature 75 mm scale length, $\pm 2.5\%$ accuracy, 5-way binding posts, and external zero adjustment
It is a laboratory apparatus for measurement of ac current in electric circuit for use in physics lab experiments
It is a analog desk type meter for the range 0-5 amperes with minimum error
Range: 0-5A.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003