

Product Name :
AC/DC Portable Suction Machine EcoVac 1

Product Code :
NLABMEDICAL-M3120001



Description :

AC/DC Portable Suction Machine EcoVac 1

Technical Specification :

Portable Suction Machine EcoVac is made by ABS Plastic with Stainless Steel back

Portable Suction Machine EcoVac is fitted with Reuseable Filter Jars : 1 x 500 ml.

The Polycarbonate Jar is fitted with lids on the side of the Suction Unit with over flow protection mechanical type

Tubing is non collapsible

Vacuum Gauge has 50 mm diameter graduated from 0 – 760 mmHg.

Bacterial Filter are Autoclavable and Reusable

Pump is Oil Free Rocker Piston Pump

Vacuum Capacity is -650 mm Hg \pm 10% at 25-30 LPM which is controlled by knob.

Noise Level<50 dB A>

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003