

Product Name :
Accelerated Curing Tank

Product Code :
NLAB-CONSTLAB270002



Description :

Accelerated Curing Tank

Technical Specification :

Driven Type: Electric

Display: Digital

Material: Stainless Steel

Operating Mode: Automatic

Size: 70.6 to 150 mm

Temperature Range (Degree Celsius): 35 deg C or 100 deg C +/- 2 deg C / 27 deg C

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

