

Product Name : Acceleration of a Geared System	Product Code : NLAB-TECHNICALAB35010
Description :	
Acceleration of a Geared System	
Technical Specification :	
Floor standing, to investigate the inertia and acceleration characteristics of a geared system. Four parallel shafts; each fitted with a sensor connected to microprocessor control unit, set of gears and inertia discs. Dead weight driven system. Safety cabinet. Set of weights, software to monitor and process the collected data. Power supply 230V/50Hz.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003