

Product Name :
Accessories to be used in conjunction with Hydraulic Bench

Product Code :
NLAB-ENGINEERINGLB19002



Description :

Accessories to be used in conjunction with Hydraulic Bench

Technical Specification :

Accessories to be used in conjunction with Hydraulic Bench Lab Equipments Manufacturer, Suppliers & Exporters

1] Eighteen accessories: Dead weight pressure gauge calibrator
Hydrostatic pressure Flow over weirs utilising bench top channel
Metacentric height Bernoulli's theorem demonstration Impact of a jet
Orifice and free jet flow Orifice discharge Energy losses in pipes
Flow channel with models Osborne Reynolds' demonstration
Flow meter demonstration Energy losses in bends and fittings
Free and forced vortex

Hydraulic ram Demonstration pelton turbine Series/parallel pumps
Centrifugal pump characteristics

[2] When coupled to the bench, the accessories are immediately
ready for use and require no additional

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003