

**Product Name :**  
Adjustable Chair

**Product Code :**  
WC144



**Description :**

Adjustable Chair

**Technical Specification :**

Adjustable Chair  
Seat Height Adjustment from 500-800 mm  
Plastic base leg  
Revolving lab stools  
With Armrest  
Ergonomic Seat Design  
Back type: height adjustable  
Backrest is made of two pcs injection moulded frame: Backrest is made of one piece injection moulded frame  
Seat is curved.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003