

Product Name :
Adjustable Chair

Product Code :
WC144



Description :

Adjustable Chair

Technical Specification :

Adjustable Chair
Seat Height Adjustment from 500-800 mm
Plastic base leg
Revolving lab stools
With Armrest
Ergonomic Seat Design
Back type: height adjustable
Backrest is made of two pcs injection moulded frame: Backrest is made of one piece injection moulded frame
Seat is curved.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003