

Product Name :
Advance Projection Microscope

Product Code :
NLAB-SCIENTIFICEQ15001



Description :

Advance Projection Microscope

Technical Specification :

Projection Microscope is very useful for routine observation and complex type measurements of specimen magnified on screen, group observations, discussions and teaching. Equipped with a 195 mm diameter Special material graduated screen Built-in base illuminator & diaphragm which is directed to the condenser which again is reflected through a mirror to the objective and compensating eyepiece for projection on to the screen again by a reflector.

Graduated Mechanical stage is also provided on the large stage for backward and forward movements.

within 30mm & 40mm scales read by verniers to 0.1mm Screen-head is rotatable through 360° with cross scales vertically and horizontally within 140mm A quadruple nose-piece carries has four objectives : 5X, 10X, 20X & 40X giving magnification from 125X to 1000X.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments