

**Product Name :**  
Advanced Infinity Corrected Binocular Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730109

**Description :**

Advanced Infinity Corrected Binocular Microscope

**Technical Specification :**

Revolving reversed quintuple objective nosepiece accepts five objectives: 4X, 10X, 20X, 40X, 100X.

Built-in double layer mechanical stage; Smooth X/Y calibrated movement; Graphite surface, rounded edges

Cord holder for convenient storage.

Equipped with a 5W LED lamp with light intensity control and a standard Kohler illumination system provides uniformly bright and free-from-glare specimen illumination.

0.1mm per fine rotation, fine focus graduation of 1um; Low position coaxial coarse and fine focus knobs.

Ergonomically designed metal frame for stability and durability. Carrying handle integrated in the back of the stand for easy carrying. Arm and base made of one piece construction.

High eyepoint, with removable eyeguards enables easy observation with/without wearing glasses.

Antifungal treatment prevents fungal growth.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003