

<b>Product Name :</b> Advanced Laboratory Metallurgical Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ173056
<b>Description :</b> Advanced Laboratory Metallurgical Microscope	
<b>Technical Specification :</b> Viewing Head: Ergonomic Sidentopf Binocular/Trinocular Head  Mechanical Stage: Right/left hand control, hard coated surface, double layer, low position coaxial movement  Body: Rugged, single die cast, ergonomic design, touch point treated  Optical system: Infinite optical system (IOS), Color corrected (CCIS), anti-fungal, anti-bacterial coated (multi-layered hard coated)  Condenser: Abbe with iris diaphragm & filter slot.  Eyepiece: Extra wide field, high point with diopter adjustment  Koehler's Illumination: Koehler halogen/LED, intensity adjustable, battery backup, auto cut-off  Incident Illumination: LED/halogen with aperture and field diaphragm  Objectives: Infinity Semi Plan/Plan/Flotar objectives  Focusing: Feather Touch Ultra Precise Coaxial (low position), both side, Coarse & Fine with tension control & upper limit stopper  Motorized: Z-axis motorized focussing  Power Supply: 90V-240V/50Hz AC  Microscopy Techniques: Polarizing, darkfield  Nosepiece: Quadruple/Quintuple with ball-bearing, slot for polarizer & analyser, rubber/metal grip & click stop	

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003