

| | |
|---|---|
| Product Name : Advanced Microscopes | Product Code : NLAB-EDUSCHOOLEQ173042 |
| | |
| Description : | |
| Advanced Microscopes | |
| Technical Specification : | |
| Features: | |
| Available as binocular and trinocular versions | |
| Nosepiece with four objectives | |
| Excellent and robust laboratory microscopes for all common routine applications | |
| Fix pre-centred Koehler illumination; adjustable field diaphragm; height adjustable Abbe condenser with adjustable aperture diaphragm | |
| A variety of optional accessories | |
| Large stage | |
| 3W dimmable LED illumination | |
| Wide field eyepieces; with achromatic infinity corrected E-plan objectives | |
| Complete with dioptre adjustment on one tube | |
| Specifications: | |
| Objective type: Infinity E-Plan | |
| Eyepieces: WF 10x / Ø 20 mm | |
| Objectives: 4x/10x/40x/100x | |

Illumination: Dimmable 3W LED (transmitted)

Tube: Binocular, Siedentopf, 30° inclined, 360° rotatable

Tube: Trinocular, Siedentopf, 30° inclined, 360° rotatable

Nosepiece: Quadruple

Interpupillary distance: 50 – 75 mm

Condenser: Abbe N.A. 1,25 pre-centred

Net weight approx.: 6.5 kg

Overall dimensions (W x D x H): 395 x 200 x 380 mm

Dioptre adjustment: Yes, one-sided

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003