

**Product Name :**  
Advanced Microscopes

**Product Code :**  
NLAB-EDUSCHOOLEQ173042

**Description :**

Advanced Microscopes

**Technical Specification :**

Features:

Available as binocular and trinocular versions

Nosepiece with four objectives

Excellent and robust laboratory microscopes for all common routine applications

Fix pre-centred Koehler illumination; adjustable field diaphragm; height adjustable Abbe condenser with adjustable aperture diaphragm

A variety of optional accessories

Large stage

3W dimmable LED illumination

Wide field eyepieces; with achromatic infinity corrected E-plan objectives

Complete with dioptre adjustment on one tube

Specifications:

Objective type: Infinity E-Plan

Eyepieces: WF 10x / Ø 20 mm

Objectives: 4x/10x/40x/100x

Illumination: Dimmable 3W LED (transmitted)

Tube: Binocular, Siedentopf, 30° inclined, 360° rotatable

Tube: Trinocular, Siedentopf, 30° inclined, 360° rotatable

Nosepiece: Quadruple

Interpupillary distance: 50 – 75 mm

Condenser: Abbe N.A. 1,25 pre-centred

Net weight approx.: 6.5 kg

Overall dimensions (W x D x H): 395 x 200 x 380 mm

Dioptre adjustment: Yes, one-sided

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003