

Product Name :
Advanced Spectrometer (Double Vernier)

Product Code :
LSP-IN-0002



Description :

Advanced Spectrometer (Double Vernier)

Technical Specification :

- **Scale:** The 177mm diameter circle is fixed, and the table and telescope movements are completely independent of each other.
- **Collimator:** Mounted on fixed pillar and fitted with 178mm focus achromatic objective, 32mm clear aperture.
- **Telescope:** Mounted on moveable pillar and fitted with 178mm focus achromatic objective, 32mm clear aperture.
- **Prism Table:** The table is marked with lines to assist placing the prism with respect to the table leveling screws, and has interchangeable clamping units for the prism and a diffraction grating.
- **Standard Accessories Supplied:** 1 Dense Flint Glass prism, 1 Prism clamp, 1 Diffraction grating holder, 1 Watchmaker's eyeglass , 1 Tommy bar / L Key for adjusting optical axis of tubes, 1 Polished hardwood case.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773