

**Product Name :**  
Advanced Student and Clinical Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ173064

**Description :**

Advanced Student and Clinical Microscope

**Technical Specification :**

Nosepiece: Quadruple/Quintuple with ball-bearing, rubber/metal grip & click stop

Viewing Head: Monocular/Binocular/Trinocular

Eyepiece: Extra wide field with diopter adjustment

Body: Rugged, Single die cast, ergonomic design, touch point treated

Koehler's Illumination: Fixed/Koehler halogen/LED, intensity adjustable, battery backup, auto cut-off, illumination mirror

Focusing: Coaxial (low position), both side, Coarse & Fine with tension control & upper limit stopper

Objectives: DIN Semi Plan/PLAN Infinite Achromatic

Power Supply: 90V-240V/50Hz AC SMPS/transformer inbuilt voltage stabilization

Mechanical Stage: Rackless, refocusing, right/left hand control, double layer, single/double slide holder, hard coated surface, low position coaxial movement

Condenser: Abbe with iris diaphragm & filter slot. Brightfield, darkfield, phase.

Optical system: Infinite optical system (IOS), colour corrected, anti-fungal, anti-bacterial coated (multi-layered hard coated)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

