

Product Name : Air Compressor Ventilators	Product Code : NLABMEDICAL-M3390007
Description :	
Air Compressor Ventilators	
Technical Specification :	
<p>A compact medical grade air compressor intended for the delivery of clean and dry compressed air.</p> <p>Adopting highly efficient oil-free compressor core and multi-level filtration and dehydration design.</p> <p>Our compressors are compatible with all the Ventilators (with external Air Source) and are very inexpensive.</p> <p>Our compressors have been running successfully with all the leading brands of Ventilators like Macquet.</p> <p>Drager, Covedian, Viasys, Siemens, Bird, Purittan Bennett, Takaoka, SLE etc.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003