

Product Name : Air Compressor Ventilators	Product Code : NLABMEDICAL-M3390007
---	---

--

Description :

Air Compressor Ventilators

Technical Specification :

A compact medical grade air compressor intended for the delivery of clean and dry compressed air.
Adopting highly efficient oil-free compressor core and multi-level filtration and dehydration design.
Our compressors are compatible with all the Ventilators (with external Air Source) and are very inexpensive.
Our compressors have been running successfully with all the leading brands of Ventilators like Macquet.
Drager, Covedian, Viasys, Siemens, Bird, Purittan Bennett, Takaoka, SLE etc.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003