

Product Name :
Air Content Meters

Product Code :
NLAB-CEMETLAB130028



Description :

Air Content Meters

Technical Specification :

The pressure gauge is installed recessed in the head of the meter. The scale has an indication range of 0-50 volumetric %. The two models are practically identical except for the capacity. This model, 1 liter capacity, is conforming to the EN 459-2 whilst, 0.75 liter capacity is conforming to EN 413-2 standards. Both models can be supplied, on request, in the electric operated version.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003