

**Product Name :**  
Air Entrainment Meter

**Product Code :**  
NLAB -CONCRLAB120053



**Description :**

Air Entrainment Meter

**Technical Specification :**

They have a quick action clamping system. Direct pressure gauge reading to the nearest 0.1% up to 6%, 0.2% from 6 to 8 and % 0.5 from 8 to 15. It is not affected by changes in barometric pressure.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003