

Product Name :
Air Flow Tester Machines

Product Code :
NLAB-TESTINGLABM89001



Description :

Air Flow Tester

Technical Specification :

Standard Methods of Testing Flexible Cellular Materials –Slabs, Bonded, and Moulded Urethane Foam
polymeric Material, Cellular Flexible – Method of Assessment of Air Flow at Constant Pressure Drop
Test Methods for Flexible Polyurethane Foam

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003