

**Product Name :**  
Air Mattress Anti Decubitus

**Product Code :**  
NLABMEDICAL-M2590006



**Description :**

Air Mattress Anti Decubitus

**Technical Specification :**

Specification :-

Mattress is made of Pvc Mattress Anti Decubitus is also known as Anti Bedsore mattress .

It is easy in operation and is better in prevention and cure of bed sores.

It replaces the foul air & reduce growth of microbes.

Power supply: 220V/50Hz, 1A/250V,

Air output: 3.5-4 liter/min

Pressure range: 50 mmHg – 90 mmHg.

It's a perfect instrument of short-term bed ridden patients.

We are Manufacturer, Exporters and Suppliers of Air Mattress Anti Decubitus at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003