

Product Name :
Air Mattress Tube Type

Product Code :
NLABMEDICAL-M2590001



Description :

Air Mattress Tube Type

Technical Specification :

Air Bed Anti Decubitis mattresses has bed sore prevention system.

Fuse rating: FIAL250V () 230V—/50Hz; 0.06A .

Maximum Weight Carrying Capacity 120 kg .

Static Mode.

Low Pressure Indicator.

Air Bed Static Mode inflates the complete mattress which makes a comfortable sitting position when the patient has to get up.

Mattress Size : 195 L x 75 W cms, Tube Height 12.5 Cm.

We are Manufacturer, Exporters and Suppliers of Air Mattress Tube Type at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003