

Product Name :
Air Permeability Tester Machines

Product Code :
NLAB-TESTINGLABM85001



Description :

Air Permeability Tester Machines

Technical Specification :

20 cm² sample holders are supplied with the device as a standard.

5-20-25-38-50-100 cm² sample holders can be provided upon request

Possibility of decreasing pressure up to 500 pa

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003