

Product Name :
Air Porosity Tester Machines

Product Code :
NLAB-TESTINGLABM91001



Description :

Air Porosity Tester Machines

Technical Specification :

Can be controlled by the touch panel

The result of air permeability test is displayed on the touch panel, units can be shifted

Equipped with printing function

Equipped with analyzing software.

11 nozzles are included

Standard Test Heads 20cm² and 38cm²

Measuring Range 1~12000mm/s (L/m².s)

Test pressure 1 ~ 4000Pa

Max thickness of specimen ? 8mm

Measuring accuracy <+/-2%

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003