

Product Name :
Air Splint Set Of 6

Product Code :
NLABMEDICAL-M2110003



Description :

Air Splint Set Of 6

Technical Specification :

* This multifunctional nursing splint can be used for Doth the left and right limbs at the same time.

The kit including :

- * Adult arm splint
- * Adult leg splint
- * Shoulder joint splint
- * Child arm splint
- * Child leg and waist splint
- * Straps and carrying bag
- * It is made of leather.
- * The splints are suitable for X-ray and MRI.

We are Exporters and Suppliers of Air Splint Set of 6 at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003